# University of Colorado Hospital Nursing Practice Guideline Interventional Radiology Checking and Exchange of Reversal/Reaction Box

### **Related Policies and Procedures:**

Adult Code Blue Emergency Resuscitation Exercise Stress Testing Outside the Heart Center

Approved by:

Pharmacy and Therapeutics Committee

**Resuscitation Committee** Department of Radiology

Nursing Practice Guidelines Subcommittee

This policy describes the procedure for assuring standardized emergency **Description:** supplies will be consistently available and safely maintained in the Interventional Radiology Area Dept Curnty Wilhi

The manager is responsible for implementation. All clinical staff is Accountability: responsible for correct use, checking and stocking.

Reversal Reaction box is an emergency support kit to be used in Interventional Radiology for sedation reversal or contrast reactions that gues IV Con Frank Definitions:

# **Policy and Procedure:**

- A. Use of Reversal Reaction Boxes
  - 1. Reversal Reaction boxes are NOT used in place of calling a code
  - 2. Reversal Reaction boxes may be used under the direction of MDs and licensed independent practitioners.
  - 3. Emergency reversal reaction boxes must be kept closed with tamper evident seal. This seal is obtained from Pyxis or the Pharmacy.

#### B. Reversal reaction Box Contents

1. The emergency medical drugs included in the reversel reaction box are primarily intended for emergencies such as anaphylactic shock or over sedation. Reversal Reaction drugs must be kept in a locked or continuously supervised location.

2. Medications Dosage guideline cards are provided for the use of the ordering provider. The adult Medication Dosage guideline eard is reviewed by UCH pharmacy staff every 3

2. Mids in Ludel for adult patients
3. an algorithm will be kept with the box.
4 Meds in cludel. G033N12

# Checking and Exchange of Reversal-Reaction Box

### C. Supply Ordering and Delivery

1. Reversal Reaction Box General-Supplies and Pharmacy supplies for new or reactin Bre replacement boxes

a. The initial emergency tackle box and replacement boxes are supplied by Interventional Radiology

b. Pharmacy will supply the drugs for the box, by removing them from the pyxis • Fax "Tackle Box Medication Order Form" (DOD OTR11076) to the in-patient

Pharmacy Storeroom the month prior to expiration.

#### D. Reversal Reaction Box Check

### 1. Daily

a. The integrity of the lock must be checked every business day by a RN or designated staff person. If clinies are closed on weekends and holidays the lock does not have to be checked on these specific days.

b. The boxes must remain locked between checks. If the seal is broken, a RN or designated staff person must recheck the entire box to ensure that all items are present. Items in the box will be used for emergent needs only.

c. The Interventional Reversal reaction Box Checklist form will be used to document daily checking of the emergency supplies box. The number on the lock and the number on the Checklist MUST be the same. If a new lock is placed, that lock number must be indicated on the Checklist on the day it was changed.

### 2. Monthly

- a. A RN or designated staff person peing the checklist must open and check the box on a monthly basis.
- b. The RN or designated staff who conducts the monthly Reversal Reaction box check will check the expiration date on the drugs. Expired drugs will be returned to pharmacy and replaced with new ones.
- c. If an item is missing, the RN or designated staff person will obtain the item from the Pyxis machine or Pharmacy via the tackle box Medication Order form and replace it.
- d. The RN or designated person will document the number of the lock placed on the emergency tackle box, on the Reversal Reaction Box checklist
- e. The Reversal Reaction Box checklist will be used to document monthly checking of the emergency box and maintained in clinical area for three years.

#### E. Post-Use

- 1. After using the box, a designated staff person will clean the box with an approved solution.
- 2. All medications used will be charged to the patient using the Pharmacy SIM number for the product(s) or Pyxis.
- 3. If charged through Pyxis, remove the replacement medication from Pyxis.
- 4. After use, replace used items from clinical supply or floor stock. Replace medications as described above. Keep the medications to restock the new box.

# Checking and Exchange of Reversal-Reaction Box

#### Table A

#### Reversal Reaction Medications

- Narcan 2mg/2ml
- Solumedrol 125mg IV
- Diphenhydramine 50mg/ml
- Elumazenil 0.5mg/5ml IV
- Atropine 1mg/10ml syringe
- Epinephrine1:10,000 1mg/10ml syringe

#### References:

- 1. Barker, L.R., Fiebach, N.H.,, Zieve, P.D., eds. *Principles of Ambulatory Medicine*, 7<sup>th</sup> Edition. Baltimore, MD: Williams & Wilkins; 2006. (LOEV).
- 2. Beers, M., Porter, R., Jones, T. et al. *Merck Manual*. 18<sup>th</sup> edition. Whitehouse Station, NJ: Merck Publishers; 2006. (LOE V).
- 3. Sommers, M., Johnson, S., Beery, T. Diseases and Disorders: A Nursing Therapeutics Manual. Philadelphia: F A Davis;2007..(LOE V)
- 4. McPhee, S., Papadakis, M., eds. *Current Medical Diagnosis and Treatment*. 49<sup>th</sup> edition. New York: McGraw Hill; 2010. (LOEV).
- 5. American Heart Association. *Guidelines for clinical exercise testing laboratories*. www.americanheart.org. Accessed August 24, 2010. (LOEV).
- 6. Auble, T., Menegazzi, J., Paris, P. Effect of out-of-hospital defibrillation by basic life support providers on cardiac arrest mortality: A metaanalysis. Annals of Emergency Medicine. 1995; 25(5):642-648. (LOE VI).
- 7. Martin, M., Lennox, P. Sedation and Analgesia in the Interventional Radiology Department. Journal of Vascular and Interventional Radiology 2003 14(9) 1119-1128 (Part I)

G033N12 Page 3 of 3